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- (1) 改訂BCOJビール分析法 8.3.6 アルコール (アルコライザー法) (2008)
- (2) AMERICAN SOCIETY OF BREWING CHEMISTS. Report of Subcommittee on Anton PAAR Alcolyzer for Measurement of Alcohol and Original Gravity. *J. Am. Soc. Brew. Chem.*, 62:195, 2004
- (3) Anton PAAR. DMA 4500/5000 Density Meter Instruction Handbook, Anton Paar GmbH, Graz, Austria, 2001.
- (4) Anton PAAR. Alcolyzer Alcohol Meter for Beer Instruction Manual, Anton Paar GmbH, Graz, Austria, 2002.
- (5) AMERICAN SOCIETY OF BREWING CHEMISTS. METHODS of ANALYSIS of the ASBC, 14th ed. Statistical Analysis-4. Youden unit block collaborative testing procedure. The Society, St. Paul, MN, 2010.
- (6) AMERICAN SOCIETY OF BREWING CHEMISTS. METHODS of ANALYSIS of the ASBC, 14th ed. Beer-4. Alcohol G. BY NEAR-INFRARED AND ORIGINAL EXTRACT CONTENT. The Society, St. Paul, MN, 2004.
- (7) Report of 2006 BCOJ Collaborative Work, Comparison of the Anton Paar alcolyzer Method and the Official GC-FID Method of the National Tax Administration Agency Japan for the Evaluation of Alcohol Content in Beer, Happo-Shu, and Nonalcoholic Beer. *J. Am. Soc. Brew. Chem.* 65, 246-247, 2007.

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- (1) AMERICAN SOCIETY OF BREWING CHEMISTS. Report of Subcommittee on Anton PAAR Alcolyzer for Measurement of Alcohol and Original Gravity. *J. Am. Soc. Brew. Chem.*, 62:195, 2004.
- (2) Anton PAAR. DMA 4500/5000 Density Meter Instruction Handbook, Anton Paar GmbH, Graz, Austria, 2001.
- (3) Anton PAAR. Alcolyzer Alcohol Meter for Beer Instruction Manual, Anton Paar GmbH, Graz, Austria, 2002.
- (4) AMERICAN SOCIETY OF BREWING CHEMISTS. METHODS of ANALYSIS of the ASBC, 9th ed. Statistical Analysis-4. Youden unit block collaborative testing procedure. The Society, St. Paul, MN, 2004.
- (5) AMERICAN SOCIETY OF BREWING CHEMISTS. METHODS of ANALYSIS of the ASBC, 2007 ed. Beer-4. Alcohol G. BY NEAR INFRARED AND ORIGINAL EXTRACT CONTENT. The Society, St. Paul, MN, 2007.
- (6) Report of 2006 BCOJ Collaborative Work, Comparison of the Anton Paar Alcolyzer Method and the Official GC-FID Method of the National Tax Administration Agency

Japan for the Evaluation of Alcohol Content in Beer, *Happo-Shu*, and Nonalcoholic Beer. J. Am. Soc. Brew. Chem. 65, 246-247, 2007.